



ANALYSEUR D'ÉLECTROLYTE ELECTROLYTE ANALYZER

Réf. : LA18070

Analyseur d'électrolyte permettant la mesure de la concentration des ions K⁺, Na⁺, Cl⁻ et Ca⁺⁺ ainsi que la valeur du pH.

- Application ISE.
- Calibration automatique en mode normal
- Mode économique.
- Imprimante thermique
- Révision automatique des données.
- Connexion RS232 et logiciel d'acquisition des données.
- Auto contrôle de l'appareil en cas d'anomalie.
- Fiable, robuste et simple d'utilisation.
- Temps de test < 1min
- Environnement de température : 10 > 30 °C
- Tension d'alimentation : AC 160V > 245V 50+ 1Hz

Electrolyte analyzer for measuring the concentration of K⁺, Na⁺, Cl⁻ and Ca⁺⁺ ions as well as the pH value.

- *ISE application.*
- *Automatic calibration in normal mode*
- *Economical mode.*
- *Thermal printer*
- *Automatic revision of data.*
- *RS232 connection and data acquisition software.*
- *Automatic control of the device in case of anomaly.*
- *Reliable, robust and easy to use.*
- *Testing time < 1min*
- *Environnement of temperature 10 > 30 °C*
- *Power supply voltage AC 160V > 245V 50+ 1Hz*

Référence	Long./Length	Larg./Width	Haut./Height	Poid/Weight
LA18070	26 cm	26 cm	34 cm	6 kg