



## ANALYSEUR DE BIOCHIMIE AUTOMATIQUE AUTO CHEMISTRY ANALYZER

Réf. : LA18040

Cet automate de chimie permet d'effectuer de manière entièrement automatique toutes les analyses courantes en chimie clinique.

- Programmation simple et intuitive (écran tactile)
- Environ 100 tests par heure
- Chimie clinique et immunoturbidimétrie
- 26 positions réactifs et 9 tubes échantillon
- 90 cuvettes
- Micro volume pour réactifs et d'échantillon
- Préchauffage réactif et détection niveaux
- Imprimante intégrée
- Automate 100% ouvert (toutes marques de réactifs acceptées).

*This chemistry analyser allows to perform automatically all the current analyse parameters in clinical chemistry.*

- *Software easy to use, simplice and flexible (touch screen)*
- *Up to 100 tests per hour*
- *Biochemistry and immunoturbidimetry*
- *26 reagents positions and 9 sample tubes*
- *90 cuvettes "on board"*
- *Micro volume for sample and reagents*
- *Reagents preheating and liquid level detection*
- *Integrated printer*
- *100% reagent open system (all reagent brands are accepted for on this device)*

Référence	Long./Length	Larg./Width	Haut./Height
LA18040	63 cm	48 cm	51 cm