

# ANALYSEUR D'HEMATOLOGIE AUTOMATIQUE 23 PARAMÈTRES 23 PARAMETERS HEMATOLOGY ANALYZER

Réf. : LA18030

**3 POPULATIONS  
23 PARAMÈTRES**  
**3 POPULATIONS  
23 PARAMÈTRES**

- Performance d'analyse : 23 paramètres, NFS formule approchée 3 populations
- Automatique
- Écran tactile

**Caractéristiques :**

- 2 canaux de comptage, 23 paramètres (incluant 3 histogrammes).
- Ecran large LCD tactile.
- 3 modes de comptage : sang complet, sang pré-dilué et sang périphérique.
- Grande mémoire, sortie multi format des résultats des tests.
- Ouverture à 3 portes pour une maintenance facile.
- Adaptateur d'alimentation externe pour réduire les bruits électroniques.

- Analytical Performance: 23 settings, NFS formula approximate 3 populations
- Automatic
- With touch screen

**Characteristics :**

- 2 counting channels, 23 parameters (including 3 histograms).
- Widescreen LCD touch screen.
- 3 counting modes: Complete blood, pre-diluted blood and peripheral blood.
- Large memory, multi-format output of test results.
- 3-door opening for easy maintenance.
- External power adapter to reduce electronic noise.

