

DIVAN D'EXAMEN EXAMINATION COUCH

Réf. : M018230

AUTRES COULEURS
SUR COMMANDE
OTHER COLORS UPON ORDER

- Structure en tube d'acier en inox rond. Relève-buste réglable. Matelas rembourré avec revêtement traité anti-feu et lavable. Porte-rouleau. Structure démontable.
- *Stainless steel tube. Fully adjustable headrest. Coating covered with washable plastic and flame retardant material. Knock down construction. Complete with paper roll holder.*



Réf.	Long./Length	Larg./Width	Haut./Height
M018230	180 cm	75 cm	55 cm

150

TABLE DE MASSAGE EN BOIS PLIABLE PORTABLE WOODEN MASSAGE COUCH

Réf. : M018240



- Table pliante en bois, hauteur réglable de 60 à 83 cm, réglage de la hauteur à différentes positions. Tiges de tension métalliques à haute résistance pour la stabilité. Revêtement en similicuir haute résistance, résistant à l'eau et à l'huile. 4 taquets en caoutchouc pour appui de la table sur le sol en position pliée. • Résistance : statique 600 kg, au travail 250 kg

- *Folding wooden massage table. Height adjustable from 60 to 83 cm / 23.4 to 32.37 in. Height adjustment in several positions. High strength metal turnbuckles for stability. High resistance eco-leather, water and oil-resistant. 4 rubber feet to support the table on the floor when folded. • Resistance: static 600 kg / 1322.8 lb. operation 250 kg / 551.2 lb*

Réf.	Long./Length	Largeur/Width	Haut./Height
M018240	182 cm	60 cm	60>83 cm