

DOUBLE POUBELLE DE BLOC DOUBLE BOWLS STAND

Réf. : BL18120



- Le porte-cupules sur roulettes est composé d'une colonne de support en tube d'acier chromé, à hauteur réglable et deux cuvettes en acier inoxydable.
- Colonne télescopique en tube d'acier chromé. Cuvette/cupule en acier inoxydable. Piètement en acier inoxydable.

• *Manufactured in chrome plated steel tube. Variable height. Stainless steel bowls.*

BAC À DÉCHETS KICK BUCKET

Réf. : BL18140



- Bac à déchets en acier inoxydable 18/10. 5 roulettes Ø50 mm. Pare-chocs.

• *Made of inox 18/10 S.S. Ø 50 mm. 5 casters with bumpers.*

CHARIOT D'ANESTHÉSIE ANAESTHESIA TROLLEY

Réf. : BL18150



- Fabriqué en acier inoxydable 18/10. Partie supérieure avec un grand espace ouvert. Six tiroirs et compartiments de rangement complets avec étagère et porte avec verrou. Monté sur roulettes diamètre 80 mm.

• *Manufactured in stainless steel 18/10. Upper section with a large open space. Six drawers and storage compartment complete with shelf and locking door. Mounted on 80 mm. casters.*

TABLE MAYO MAYO TABLE

Réf. : BL18130



- Table Mayo en acier inoxydable, hauteur variable, 4 roulettes.

• *Stainless steel construction, variable height, 4 castors*

	Long./Length	Larg./Width	Haut./Height		Long./Length	Larg./Width	Haut./Height
BL18130	450 mm	675 mm	800/1150 mm	BL18150	40 cm	60 cm	110 mm