

## POMPE À PERFUSION INFUSION PUMP

Réf. : RE18020



- Écran LCD TFT couleur de 4,3 " (rétro-éclairage LED), avec résolution 272 × 480
- Trois modes de fonctionnement : débit / temps / volume
- Bibliothèque de médicaments avec 210 types de médicaments
- Mémoire : 1500 historiques
- La fonction de chauffage, adaptée à l'infusion en hiver ou nécessite la température du médicament liquide
- Affichage multi-langue, avec 6 heures d'autonomie de batterie
- Variété d'alarme visible et sonore, ce qui rend la perfusion plus sûre
- Télécommande pour changer le mode, régler le débit et tout les paramètres
- 4.3" color TFT LCD screen (LED backlight), with resolution 272×480
- Three working modes: Rate/Time/Volume mode
- Drug library with 210 kinds of medicine
- Storage 1500 history records
- Heating function, suitable for infusing in winter or there is requirement for the liquid medication temperature
- Support multi language display, with 6 hours battery autonomy
- A variety of visible and audible alarm, make infusion more safe



## POUSSE SERINGUE "TM500" SYRINGE PUMP "TM500"

Réf. : RE18030

- Volume de la seringue: plage de 0,1 à 99,9 ml (0,1 ml)
- Volume Total : plage : 9999,9 ml (échelle de 0,1 ml), précision ± 2%
- Alarmes: de pré-fin de perfusion, fin de perfusion, volume pré-réglé, débit nul ou départ non, occlusion, débrayage mécanique, erreur de réglage, batterie faible, coupure secteur et d'autres défaillances d'appareil.
- Pression atmosphérique 86kPa ff 106kPa
- Niveau de sécurité : classe I, type BF, batterie intégrée, dispositif de fonctionnement continu
- Télécommande pour changer le mode, régler le débit et tout les paramètres
- Syringe volume: range from 0.1 to 99.9 ml (0.1 ml)
- Total Volume: 9999.9 ml (0.1 ml step), accuracy ± 2%
- Alarms: almost empty, empty syringe, volume pre-set, complete injection, occlusion, loose syringe, adjustment error, low battery, power failure and other device failures.
- Atmospheric pressure 86kPa ff 106kPa
- Safety level: Class I, type BF, built-in battery, continuous operation

Long./Length	Larg./Width
282 mm	144 mm
Haut./Height	Poids/ Weight
148 mm	2.300 kg