

UNITÉ DENTAIRE "MÉGADENTAL" "MEGADENTAL" DENTAL UNIT

Réf. : DE18010

- Fauteuil dentaire très design et ergonomique
 - Système de nettoyage du filtre dans l'aspiration de la pression négative de l'air et l'éjecteur salivaire de la pression négative de l'eau.
 - Éjecteur salivaire sous pression.
 - Système de verrouillage de pression d'air pour une manipulation facile.
 - Système moteur LINAK origine Danemark pour un mouvement stable.
 - Accoudoir côté extérieur mobile pour patient.
 - Appui-tête peut être réglé et verrouillé.
 - Couvercle anti-incrustation en plastique dans le bac à instruments, interrupteur principal d'eau et d'air
 - Alimentation : AC 220V. 50H z / 60H z
 - Alimentation de l'air : pression de l'eau 200-400 kPa
 - Alimentation de l'air : pression atmosphérique 550-800 kPa
- *Very design and ergonomic dental chair*
 - *Filter cleaning system installed in Air negative pressure suction and water negative pressure salivary ejector.*
 - *Air pressure locking system for easy manipulation.*
 - *Adopt LINAK motor system from Denmark for stable movement.*
 - *The out side armrest for patientis mobile.*
 - *The headrest could adjust and lock.*
 - *Plastic Antifouling Cover the in instrument tray, main water and air switch.*
 - *Power Supply: AC 220V. 50Hz/60Hz*
 - *Air source: Air pressure 550-800kPa*
 - *Air source: Water pressure 200-400kPa*



Référence	Long./Length	Larg./Width	Haut./Height	N/W kg	G/W kg
DE18010	1490 mm	1020 mm	1120 mm	210	300