


## NEGATOSCOPE EXTRA PLAT EXTRA THIN X-RAY VIEWER

- La structure de l'écran est en alliage d'aluminium haute résistance. Les angles sont protégés par une plaque d'acier, recouverte d'ABS. Ecran en polymère à méthacrylique. Le corps du négatoscope est en alliage d'aluminium haute résistance. L'intérieur et l'extérieur sont recouverts d'une poudre anti-statique.
  - Le système de fixation pour radios permet une installation et un retrait facile.
  - Voltage/Fréquence : 220V/ 50-60Hz
  - Luminosité : 3 000 cd/m<sup>2</sup>
  - Option : 3 plages
- *The screen structure is made of high-strength aluminum alloy. The corners are protected by a stainless plate, covered with ABS. Methacrylic polymer screen. The body of the viewer is made of high-strength aluminum alloy. The interior and exterior are covered with an anti-static powder.*
  - *A fastening system for the radios allows a simple installation and removal.*
  - *Voltage / Frequency: 220V / 50-60Hz*
  - *Luminance: 3.000 cd / m<sup>2</sup>*
  - *Option : 3 sections*



- 1 plage
  - 1 section
- Réf. : DA18170

- 2 plages
  - 2 sections
- Réf. : DA18171



	Long. /Length	Larg. /Width	Haut. /Height
DA18170	44 cm	2.4 cm	58 cm
DA18171	72 cm	2.4 cm	58 cm