

# DIVAN ÉLECTRIQUE “LEMI 4” COUCH ELECTRIC “LEMI 4” Réf : M018220

• Réglage totalement électrique de tous ses mouvements. Idéal pour tous les opérateurs qui veulent travailler en conditions optimales et accueillir leurs clients avec style. Parfaitement adaptable et simple à utiliser, cette table est pourvue de 4 moteurs électriques qui permettent de choisir en autonomie la position et la hauteur plus indiquées pour travailler en offrant au patient la possibilité de s’asseoir sans efforts. La centrale de commande, réalisée sous projet de Lemi Group, est pourvue d’une commande automatique de retour à zéro, qui ramène le fauteuil à sa position initiale

**Norme :**

- 4 moteurs :  
(hauteur, dossier, repose-jambes, trendelenburg)
- Contrôle du combiné

• All-electrical adjustment of all its positions, Perfect companion for any operator wanting to accommodate his/her customers in style, offering excellent working conditions. Extremely versatile and easy to use, the bed comes complete with 4 electric actuators which allow the operator to autonomously select the most suitable working height and position, thereby providing for maximum customer accessibility. The electronic control unit, custom designed by Lemi Group, includes a useful automatic reset control which allows the bed to be brought back to its original position.

**Standard:**

- 4 motors (height, back rest, leg rest, trendelenburg)
- Handset control

| Réf.    | Longueur /Length | Largeur/Width |
|---------|------------------|---------------|
| M018220 | 162 cm           | 62 mm         |



**EN OPTION**  
**OPTIONAL**  
**Pédale**  
*Pedal*  
Réf. : M018211



**CHAISE CHAIR**  
Réf. : M018222



**TABOURET STOOL**  
Réf. : M018223



OPTIONS POUR LE LEMI 3 ET LEMI 4