

CHARIOT INOX MULTIFONCTIONS STAINLESS STEEL MULTIFUNCTIONAL TROLLEY

Réf. : M018100

• Bâti et tiroirs en acier laqué epoxy anti-moisissures. Plan de travail en technopolymères pour faciliter le nettoyage et la désinfection. Réhausse en acier inox parfaitement intégrée à la structure avec pilulier. Façade : glissières métalliques à billes silencieuses dotées de butées de fin de course, 3 tiroirs 615x470x147 mm et 1 tiroir 615x470x223 mm. Fermeture centralisée à clé. Côté droit: 1 tablette ABS escamotable, 1 poubelle 10L en ABS, 1 porte-sondes cathéters en ABS. Côté gauche: 4 espaces avec galerie fil métallique. 4 roues pivotantes Ø125mm, anti-traces et ne s'effilochant pas, dont 2 avec frein indépendant. Galets pare-chocs PVC antitraces aux angles. Chariot en acier inox.

• *Multifunction cart with: 4 drawers in painted steel, 1 ABS extending table, 1 ABS 10 liters waste bin, 1 grey ABS catheter and probe holder. Self-bearing monocoque body made of steel sheet coated with anti-mould epoxy powders. The left compartment features four trays with metal rim. Drawers are made of anti-mould epoxy powder coating steel, and are equipped with: central locking system, silent sliding system on metal tracks and rollers, end stoppers. 3 drawers mm. 615 x 470 x 147 H, 1 drawer mm. 615 x 470 x 223 H. The tabletop is made of technopolymers and has no sharp edges for easy cleaning and disinfecting. Holding frame for dispenser with modular drawers. Push handle in stainless steel, fully integrated in the frame. 4 casters, 2 with brake, Ø 125mm, non-marking, anti-thread, mounted on swiveling bearings, 4 PVC non-marking bumpers. Stainless Steel cart.*



	Long./Length	Larg./Width	Haut./Height
M018100	62 cm	80 cm	100 cm