

TABLE DE RADIO OS POUMONS **NUMÉRIQUE** X-RAY BUCKY ROOM **DIGITAL**

Réf. : IM18100

Notre table de radiologie numérique polyvalente met à votre disposition tous les avantages liés à vos besoins, cette technologie couvre une large gamme d'examens.

- Un dosage parfaitement contrôlé grâce à un dispositif de détection et de mesure des radiations ionisantes.
- Un générateur convertisseur haute fréquence qui optimise la qualité des images et une excellente reproductibilité.
- Petite distance focale, grande capacité de masse, tube à rayons X de haute performance.
- Platine avec verrouillage électromagnétique 4 positions.
- Affichage LED numérique.



Our digital X-RAY bucky room gives you all the advantages linked to your needs and our technology. Conventional radiological table

- *A perfectly controlled dosage thanks to a detection device and of ionizing radiation.*
- *A high-frequency converter generator optimizes image quality.*
- *Small focal length, large mass capacity and high performance X-ray tube.*
- *Turntable with electromagnetic lock 4 positions.*
- *Digital LED display (light emitting diode).*

	Long./Length	Larg./Width	Haut./Height	Poids/Weight/N.g	Poid/Weight/G.g
Tableau radiographique/Radiography Table	2290 mm	1340 mm	910 mm	410 kg	546 kg
Support/Bucky Stand	2360 mm	800 mm	550 mm	140 kg	240 kg
Console de commande/Control console	1760 mm	1085 mm	915 mm	110 kg	133 kg
Générateur HV/HV Generator	986 mm	547 mm	1077 mm	120 kg	146 kg

TABLE DE RADIO OS POUMONS **ANALOGIQUE** X-RAY BUCKY ROOM **ANALOGIC**

Réf. : IM18101

