

TABLE À LANGER CHANGING TABLE

Réf. : MA18140

- Structure en tube d'acier, revêtement époxy. Plan de couchage revêtu de vinyl lavable. Une étagère. Structure démontable.

- *Made of epoxy painted steel tube. Top and sides upholstered and covered with washable plastic material. Laminated plastic lower shelf. Knock down construction.*



Réf.	Long. /Length	Larg. /Width	Haut. /Height	Capacité/ Capacity
MA18140	177 cm	66 cm	90 cm	50 kg

66

VENTOUSE OBSTETRICALE MECANIQUE MECHANICAL OBSTETRIC MACHINE

Réf. : MA18150



- **Comprenant :**
3 ventouses complète diam. : 40, 50 et 60 mm, 1 support flacon en fil plastifié, 1 flacon d'aspiration en verre, 1 pompe manuelle pour dépression, 1 bouchon avec manomètre 1 jeu de tuyaux caoutchouc

- **Including:**
3 complete suction cup diameter : 40, 50 and 60 mm, 1 plastic bottle holder 1 glass suction bottle, 1 manual vacuum pump 1 cap with manometer, 1 set of rubber hoses

PÈSE BÉBÉ BABY SCALE

Réf. : MA18160

Capacité

- Poids minimum 0,5 kg
- Poids maximum 50 kg
- Graduation 0,01 kg
- Unité de poids kg ou lb-oz
- Tolérance 5‰
- Alimentation 2 piles AA (fournies)

Capacity

- *Minimum weight 0.5 kg*
- *Maximum weight 50 kg*
- *Graduation 0.01 kg*
- *Weight unit kg or lb-oz*
- *Tolerance 5 ‰*
- *Power 2 AA batteries (included)*

