

## AUTOCLAVE 365 LITRES AUTOCLAVE 365 LITERS

Réf. : ST18010

**365  
LITRES  
LITERS**



### Caractéristiques :

- Fabriqué en acier inoxydable AISI 304
- Fermeture automatique et système de verrouillage opéré par un cylindre pneumatique
- Mouvements vertical et descendant
- Isolation thermique
- Scellement à travers l'injection d'air comprimé et la garniture de silicone
- Capteur photoélectrique de sécurité anti-écrasement
- Capteur d'interverrouillage de porte qui empêche l'ouverture simultanée des deux portes, évitant ainsi la contamination croisée

### Features:

- Made of AISI 304 stainless steel or higher
- Automatic closing and locking system operated by a pneumatic cylinder
- Vertical ascendant and descendant movement;
- It has thermal insulation
- Sealing through compressed air injection, through siliconetrims
- Safety anti-crushing photoelectric sensor
- Door interlocking sensor that prevents the two doors to be opened together, thus avoiding cross contamination

Réf.	Dimensions interne/Internal dimensions			Dimensions externe/External dimensions		
	Long. /Length	Larg. /Width	Haut. /Height	Long. /Length	Larg. /Width	Haut. /Height
ST18010	66 cm	46 cm	120 cm	204 cm	94 cm	169 cm