

MALLETTE MÉDICAL D'URGENCE EMERGENCY MEDICAL CASE

Réf. : UR18080

- Mallette d'urgence en aluminium. La mallette peut contenir différents types d'équipements médicaux. En option : aspirateur RES-Q-VA
- *Emergency medical case is ideal for doctor's first aid and emergency operation, made of anodized aluminium, it is resistant and ensure long durability. Optional: suction pump RES-Q-VAC.*

Réf.	Long. /Length	Larg. /Width	Haut. /Height
UR18080	52 cm	34 cm	16 cm



BRANCARD PLIABLE FOLDABLE STRETCHER

Réf. : UR18090

- Brancard d'urgence pliable, bâti en aluminium anodisé, toile de polyester imperméable, lavable. Pieds en aluminium. Capacité 120 kg.
- *Foldable emergency stretcher, made of anticorrosive and anodised aluminium. With a polyester, fire retardant and anty putridity sheet. Aluminium feet. Capacity 120 kg.*



Réf.	Long. /Length	Larg. /Width	Haut. /Height	Poids /Weight	Surface utile /Useful surface
UR18090	220 cm	52 cm	12 cm	7 kg	173 x 52 cm

RAPPEL REMINDER

PAGE 123

OXYMÈTRE DE POULS DIGITAL PULSE OXYMETER

Réf. : DA1816 0

